



DR. **MÜLLER**  
FOREVER YOUNG



# DR. MÜLLER BEAUTY LIGHT THERAPY 64

Product Leaflet 2020-2021

# DR. MÜLLER BEAUTY LIGHT THERAPY 64

## LAMP SPECIFICATIONS

Total lamps

24 x Dr. Müller Collagenic Light 200W E-Tronic | 200 cm  
 20 x Dr. Müller Bright Skin Light 200W E-Tronic | 200 cm  
 20 x Dr. Müller Blue Skin Light 200W E-Tronic | 200 cm

Recommended exposure time

10 minutes

## TECHNICAL SPECIFICATIONS

Power supply

3 x 400V + N + PE

Plug and cable

32A CEE 5 x 4 mm<sup>2</sup>

Power

12.400W

Breaker/fuse (slow)

3 x 25A

Minimum net weight<sup>(1)</sup>

290 kg

Minimum crated weight<sup>(1)</sup>

335 kg

Dimensions open

176 x 162 x 227 cm

Dimensions closed

162 x 151 x 227 cm

Boxed dimensions

227 x 114 x 156 cm

## + AIRCO

3 x 400V + N + PE

32A CEE 5 x 6 mm<sup>2</sup>

13.400W

3 x 35A

350 kg

335 kg + 90 kg

210 x 189 x 238 cm

199 x 162 x 238 cm

227 x 114 x 156 cm

+ 120 x 80 x 90 cm

## COLOUR OPTIONS



Dream White Metallic



Fancy Red Metallic

## STANDARD EQUIPPED WITH



IQ Touch Control



Age Control



Electronic Ballasts



TopCooler



Rainbow Manager



VibraPlate



Music Pro



Wellness Sound



Voice Guide



Bluetooth



Handles



Air Extraction Hose

## ADDITIONAL



Breeze



Aroma



Airco



Changing Cubicle

## CONTROL SYSTEMS



Dr. Müller Manager



Paytronic



Wireless Management System

<sup>(1)</sup> Minimum weight is based on a Dr. Müller Beauty Light Therapy 64 with standard configuration. Machine specific weight information available upon request.

When placing the Dr. Müller Beauty Light Therapy 64 make sure that the device is surrounded by 10 -15 centimetres of free space for air circulation.



*IQ Touch Control*



*Music Pro*



*TopCooler*



# FOREVER YOUNG



Granaatstraat 6 | 7554 TR Hengelo | The Netherlands  
+31 (0)74 2912215 | [info@drmuller.com](mailto:info@drmuller.com) | [DrMuller.com](http://DrMuller.com)